

E-Z HEAT WARMING SYSTEM

E-Z Heat is an extremely user-friendly system that allows up to four connected beds to help maintain the patients' body temperatures during surgical procedures. This patented warming system allows the user to heat the induction chamber and warming beds simultaneously. The warming beds range from 95 to 104 degrees Fahrenheit, or 35 to 40 degrees Celsius. Bed temperatures can be adjusted in the administrator options. Each of the warming circuits displays the elapsed time and the set temperature, which can be individually controlled. All settings made in the administrator options are password-protected to ensure security from tampering individuals. The heating bed also allows the option to mount our Microflex Breather and works in perfect harmony with all E-Z Anesthesia systems.



E-Z Heat System Features:

- 7" Touch screen control display
- Stainless steel heating beds 10" L x 4" W
- Multiple languages available
- Very user-friendly
- Remote access for added support
- Voltage input: 120-240 VAC/50-60 Hz



EZ-HEAT Controller



HB-163XL Heating Bed

Use as chamber warmer or surgical bed for larger animals.

Surgical Heating Beds to Suit Your Needs

Choose from Flat, V-shape, or XL Flat beds. Warming beds are constructed of heavy gauge 304 stainless steel. Heat evenly distributes throughout the bed surface. Beds range from 95 to 104 degrees Fahrenheit or 35 to 40 degrees Celsius.



HB-163XL Heating Bed
(11" L x 6" W x 1" H)

HB-1000-V V-Shaped Surgical Heating Bed
(10" L x 4" W x 1" H)

HB-1000 Flat Surgical Heating Bed
(10" L x 4" W x 1" H)

Add Breathers

Optional EZ-103H Microflex Breather easily attaches to bed with bracket.



HB-103XL Heating Bed with Breather
(11" L x 6" W x 1" H)

HB-103V V-Shaped Heating Bed with Breather
(10" L x 4" W x 1" H)

HB-103 Flat Heating Bed with Breather
(10" L x 4" W x 1" H)